

# BIODERM ESSENTIAL : VACUUM ESSENTIAL OIL SYSTEM

**VACUUM THERAPY ESSENTIAL OIL SYSTEM** combines two well-proven techniques : **Essential Oil System (patented system)** and **mechanical palpate-roll**. The combination of those two techniques, in a global slimming concept, offers good results particularly in cases where excess weight is hard to move, and orange peel problem.

As a combination of the mechanical effects of vacuum therapy and the chemical properties of essential oils, the **BIODERM Essential** enables efficient, comfortable and unique treatments.

**The technique :** the new patented system can manifest itself through the spraying of essential oils through the roller head during the "palpate-roll". Absorbed by the skin, oils penetrate the circulation (osmotic power), they stimulate the nervous system by a modification of the chemical structure of body liquids and tissue. The essential oils favor the body self-regulation. Their properties are various: analgesic, toning, disinfectant, anti-inflammatory, draining... They provide well being feelings, especially thanks to their olfactory property. The slide of the head on the skin makes actions effectuated on the skin comfortable and reinforces the efficiency thanks to the combination of agents. The beauty and the suppleness of the skin at the end of the treatments, the perfumes exhaling from the oils in the cabin participate to the customer well being.

In the area of aesthetic treatments, the physiological effects encouraged refer to the characteristics of the dermis and the hypodermis:  
Vascularization • Defibrosing • Drainage • Reflex effects • Relaxation

**Face treatment :** Thanks to the use of the proper massage head, the Bioderm technique can be applied to the face; this is called Ridoplastie. Wrinkles are the result of skin ageing and also fibrinolytic phenomena which attach the base of the wrinkle to the tissues.

The face and neck are areas with lots of lymphatic and venous exchanges and are often insufficiently stimulated. Ridoplastie face cares give noticeable results from the first session as regards both external aspect and feeling of well-being.

**BIODERM Essential** enables to put in place treatments as well as for the scars or the labor strike.



THE METHOD

BODY ASSESSMENT

NUTRITIONAL ASSESSMENT

TREATMENT  
OF FIBROSE TISSUES

VASCULARISATION  
DEFIBROSING  
DRAINAGE

LYMPHATIC DRAINAGE

ELECTROSTIMULATION



**The equipment :**

- Bioderm offers free or programmed work methods; it can be used in continuous or pulsated mode; the pressure goes from 50 to 900 mbar.
  - 3 body draining programmes: trunk-belly, legs-thighs, general back.
  - 3 slimming programmes: waist-stomach, thighs, ankles.
  - 4 face-cares programmes: face top and bottom, lifting, wrinkles.
  - 3 scars treatment programmes: face top and bottom, body.

- Continuous mode – pulsated mode: Frequency: 0,5 Hz/0,7 Hz/1 Hz/2,5 Hz.
- Mains voltage: 220/240 V to 50/60 Hz.

Contra-indications : skin tumours, inflammatory rheumatism, major dermatosis, progressive infectious diseases, tumours, excessive hair fragility.



**Massage heads**

- Head 1 & 2: Scars and wrinkles treatment
- Head 3: scars and body treatment
- Head 4-5-6: body treatment.

