

SP3+[®] ELECTROMYOTHERAPY

The latest electro-therapy device of Multiwell enables two techniques of work: Lipo-stimulation and Myo-stimulation

- Lipo stimulation technique enables to treat the cellulite in releasing the fat acid by sending a painless current through the hypodermis.
- Myo-stimulation offers a passive muscular activity by sending a painless current through the muscle.

Its 12 outputs enable to work with 24 electrodes, to treat several areas in the same time and to combine lipo and myo-stimulation program together.

The type of current used, its form, its intensity, its frequency (and other parameters) must be adapted to the muscles chosen and the aims: body toning and remodeling, strength or volume increase, vascularisation, relaxation...

SP3+ can also be used as a complement for physical preparation, muscular training for injured people, and also simply for persons who do not want to practice sports.

The power of the SP3+ will give to your patients a full reshaping and an efficient muscular toning up without any effort.

- **Performances and time saving** : during the same session, thanks to the 24 electrodes, it is possible carry out different programmes on different parts of the body; and therefore both slimming treatment on an area of excess fat and the same time as muscle toning (electrostimulation) on one or two muscle groups. For instance, it is possible to work simultaneously on the hips, the buttock and the abdominal muscles.

- **Comfort** : SP3 is equipped with a unique concept whereby the current installation gradient can be adjusted so that the different muscle fibres are gradually called on making stimulation more comfortable.

- **Easy to use** : SP3 is equipped with specific programmes for each part of the body or each type of use. Just select the right programme and adjust the intensity level to that required for the treatment:

- Specific slimming programmes
- Mixed slimming and muscle toning programmes
- Specific programmes for each muscle group: buttocks, abdominal muscles, pectoral muscles; back muscles, quadriceps, etc...
- Active recovery programmes for volume, explosive strength, resistance, endurance, warming-up, relaxation, vascularisation.
- Specific programmes for face treatment.
- 600 possibilities of programmation
- Easy Start programming



Safety :

SP3 generators are equipped with secured outlets to protect the customer from any risk of electric shock or burning.

Characteristics :

Biphasic current of nil average.
Power supply:
220V/50Hz (alternative sinusoidal)
Option : 110 V
Absorbed power :
65W Device of type IIBF
Fuses : 2 fuses 240V/0.3A 4x20mm fast fusion
In conformity with directive 93/42/CEE



THE METHOD

BODY ASSESSMENT

NUTRITIONAL ASSESSMENT

TREATMENT
OF FIBROSE TISSUES

VASCULARISATION
DEFIBROSING
DRAINAGE

LYMPHATIC DRAINAGE

ELECTROSTIMULATION